

Walking frames and Rollators



www.rebotec.de

Walk-on

Aluminium walking frame

Equipment:

- powder coated aluminium frame
- nonslipping handles for safer use
- foldable, height adjustable
- approved up to 150 kg body weight

ref. no. 174. ...

Up to 150 kg!



ref. no. 174.10



ref. no. 174.40



Type:	XS	S	M	L
depth:	43 cm	43 cm	52 cm	52 cm
width:	52 cm	52 cm	59 cm	59 cm
height of handle:	62 cm - 71 cm	70 cm - 81 cm	68 cm - 87 cm	78 cm - 97 cm
weight:	3,0 kg	3,1 kg	3,6 kg	3,7 kg
adjustments in height:	5x	5x	10x	10x
max. body weight:	150 kg	150 kg	150 kg	150 kg
colour:	coloured	coloured	anthracite	anthracite
ref. no.:	174.10	174.20	174.30	174.40

Walk-on

Aluminium walking frame with 2 rolls



ref. no. 175.40

Equipment:

- powder coated aluminium frame
- nonslipping handles for safer use
- foldable, height adjustable
- approved up to 120 kg body weight

ref. no. 175. ..



ref. no. 175.10

Type:	XS	S	M	L
depth:	40 cm	40 cm	47 cm	47 cm
width:	50 cm	50 cm	59 cm	59 cm
height of handle:	62 cm - 71 cm	70 cm - 81 cm	68 cm - 87 cm	78 cm - 97 cm
weight:	3,3 kg	3,4 kg	4,0 kg	4,1 kg
adjustments in height:	5x	5x	10x	10x
max. body weight:	120 kg	120 kg	120 kg	120 kg
colour:	coloured	coloured	anthracite	anthracite
ref. no.:	175.10	175.20	175.30	175.40

Fixi and Fox

Rollator for children

Equipment:

- powder coated aluminium frame
- dual action brake
- foldable
- airless tyres
- approved up to 100 kg body weight

re. no. 186.00. ..

Accessories:

- seat (black)
ref. no. 286.65.25
- back pad
ref. no. 280.10.48

Only 6,2 kg!



ref. no. 186.00.10

Type:	Fixi	Fox
depth:	68 cm	68 cm
width:	58 cm	58 cm
height of handle:	61 cm - 70 cm	70 cm - 78 cm
seat height (optional):	47 cm	47 cm
seat width (optional):	38 cm	38 cm
weight:	6,2 kg	6,3 kg
adjustments in height:	5x	5x
max. body weight:	100 kg	100 kg
ref. no.:	186.00.10	186.00.20



Pluto

Children reverse walker



ref. no. 187.00.10
shown with seat size 1

Equipment:


- 2 different sizes
- width and height adjustable handgrips
- return motion break
- steerable frontwheels
- airless wheels

ref. no. 187.00. ...

- the Pluto reverse walker provides independent movement for children
- it supports the user from behind and promotes upright posture when standing and walking
- it supports the therapeutic exercises to learn natural movements and to increase the balance
- the braking system in the back wheels prevents unintentional rolling back

Accessories:

- front-wheels with steerlock
ref. no. 287.10.42
- seat for Pluto - size 1: up to 60 kg
ref. no. 287.20.10
- seat for Pluto - size 2: up to 80 kg
ref. no. 287.20.20



Type:	Size 1	Size 2
length:	79 cm	85 cm
width:	64 cm	67 cm
height of handle:	45 cm - 81 cm	55 cm - 91 cm
weight:	8,0 kg	8,5 kg
max. body weight:	60 kg	80 kg
ref. no.:	187.00.10	187.00.20

Polo Plus-T

Aluminium rollator with bag

Equipment:

- powder coated aluminium frame
- dual action brake
- foldable
- seat reversible to tray
- incl. bag
- airless tyres
- approved up to 120 kg body weight

ref. no. 181...

Accessories:

- back pad
ref. no. 280.10.48
- holder for walking stick
Art.-Nr. 164.30.00
- removable shopping bag
ref. no. 284.00.70



ref. no. 181.20

*Available in
3 seat heights!*

Type:	S	M	L
depth:	70 cm	70 cm	70 cm
width:	67 cm	67 cm	67 cm
height of handle:	64 cm - 74 cm	69 cm - 82 cm	74 cm - 92 cm
seat height:	49 cm	55 cm	60 cm
seat width:	48 cm	48 cm	48 cm
weight:	7,4 kg	7,6 kg	7,8 kg
adjustments in height:	5x	7x	10x
max. body weight:	120 kg	120 kg	120 kg
ref. no.:	181.10	181.20	181.30

Polo-Comfort

Rollator fully equipped



ref. no. 182.20

*Available in
2 seat heights!*



Separately removable
shopping bag.

Equipment:

- powder coated aluminium frame
- dual action brake
- back pad
- foldable
- seat reversible to tray
- incl. bag with removable shopping bag
- holder for walking stick
- airless tyres
- approved up to 120 kg body weight

ref. no. 182. ...



Type:	M	L
depth:	70 cm	70 cm
width:	67 cm	67 cm
height of handle:	69 cm - 82 cm	74 cm - 92 cm
seat height:	55 cm	60 cm
seat width:	48 cm	48 cm
weight:	7,8 kg	8,0 kg
adjustments in height:	7x	10x
max. body weight:	120 kg	120 kg
ref. no.:	182.20	182.30

Jumbo 150

Steel rollator



Equipment:

- powder coated steel frame
- dual action brake
- foldable
- seat reversible
- incl. bag
- airless tyres
- 10 adjustments in height
- approved up to 150 kg body weight

ref. no. 188.00.00

Accessories:

- back pad
ref. no. 280.10.48
- holder for walking stick
ref. no. 164.30.00



ref. no. 188.00.00
shown with back pad

Type:	Jumbo 150
depth:	72 cm
width:	76 cm
height of handle:	75 cm - 93 cm
seat height:	60 cm
seat width:	54 cm
weight:	11,6 kg
adjustments in height:	10x
max. body weight:	150 kg
ref. no.:	188.00.00

Up to 150 kg!